



exotech[®]

MILK ADULTERATION TEST KIT

Milk is a nutrient-rich beverage that is packed with important vitamins and minerals, including calcium, phosphorus, B vitamins, potassium, and vitamin D. It is also an excellent source of protein. Overall, milk is a healthy beverage that can provide a number of health benefits. If you are able to tolerate milk, it is a good choice to include in your diet.

HUMAN HEALTH

Milk adulteration is the practice of adding harmful or substandard substances to milk in order to increase its volume or improve its appearance. The health effects of consuming adulterated milk can vary depending on the type of adulterant that is used.


EXOTECH MILK ADULTERATION TEST KIT

Exotech is a revolutionary solution to combat milk adulteration, providing a practical and effective means to ensure the purity and quality of milk. By empowering individuals with the ability to detect adulteration, Exotech aims to restore trust in the dairy industry and protect the health of consumers. With its user-friendly design, accurate results, and affordability, Exotech is set to revolutionize the way milk is tested and consumed, creating a safer and more transparent dairy market.

exons.pk@gmail.com
www.exonscientific.pk




STARCH DETECTION TEST




Procedure	Add 2-3 drops of reagent in 3 ml milk hot milk	
Detection Limit	500 ppm (0.05 %)	
Control		Creamy/Light Yellow
Positive		Light Blue with Black Particles

HYDROGEN PEROXIDE DETECTION TEST



Procedure	Add 0.3 ml of reagent in 2 ml milk	
Detection Limit	750 ppm (0.075 %)	
Control		Light Grey
Positive		Light Blue to Greyish Blue

NEUTRALIZERS DETECTION TEST



Procedure	Add 1 ml of reagent in 1 ml milk	
Detection Limit	150 ppm (0.015 %)	
Control		Light Orange
Positive		Light Pink to Dark Pink

BORIC ACID DETECTION TEST




Procedure	Add 1 ml R1 and 1 ml R2 reagents in 1 ml milk	
Detection Limit	1000 ppm (0.15 %)	
Control		Dark Pink Color
Positive		Appearance of White

COMMON SALT DETECTION TEST




Procedure	Add 1 ml R1 and 0.25 ml R2 reagents in 2.5 ml milk	
Detection Limit	750 ppm (0.075 %)	
Control		Dark Brown
Positive		Distinct Yellow

SMP /WHEY POWDER TEST




Procedure	Add 0.25 ml R1, 1 ml R2, 1.5 ml milk and put in boiling water bath for 5 min	
Detection Limit	30,000 ppm (3.0 %)	
Control		Light Green
Positive		Greyish to Bluish Grey

MILK UREA DETECTION TEST



Procedure	Add 1 ml of reagent in 1 ml milk	
Detection Limit	750 ppm (0.075 %)	
Control		Light Yellow
Positive		Distinct Yellow

FORMALIN DETECTION TEST *



Procedure	Add 1 ml of reagent in 1 ml milk	
Detection Limit	25 ppm (0.0025 %)	
Control		Light Brown
Positive		Purple Ring

DETERGENT DETECTION TEST




Procedure	Add 500µl R1 and 1ml R2 in 500µl milk and centrifuge for 3min at 3000 RPM	
Detection Limit	750 ppm (0.075 %)	
Control		More Intense Blue Color in Upper Layer
Positive		More Intense Blue Color in Lower Layer

HYPOCHLORITE DETECTION TEST **



Procedure	Add 1 ml of reagent in 3 ml milk	
Detection Limit	1000 ppm (0.1 %)	
Control		Creamy White
Positive		Light Yellow to Yellow

CANE SUGAR DETECTION TEST



Procedure	Add 1 ml of reagent in 1 ml milk and put in boiling water bath for 5 minutes	
Detection Limit	750 ppm (0.075 %)	
Control		Light Brown
Positive		Red to Light Brick Red



Storage Conditions: Store at room temperature and avoid exposure of direct sunlight.

Close lid of bottle tightly after usage to prevent spillage and reagent evaporation.

* Don't shake tubes after addition of reagents.

** Read result within 5 seconds of reagent addition.